

"FINALLY A WHEELCHAIR THAT GETS ME TO THE TOP OF MY DRIVEWAY. SOMETIMES IT'S THE LITTLE THINGS LIKE BEING ABLE TO COLLECT MY OWN MAIL THAT MAKE A DIFFERENCE"



Sue
United Kingdom

magic
Mobility™

Frontier V6 Compact 40

POWER IN A COMPACT PACKAGE

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.



UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be custom designed to a fixed position. It's whatever you need.

TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.



V3/AUST/MARCH19

BUILT TO LAST

CUSTOMISED FOR YOUR LIFESTYLE



AUSTRALIAN MADE

Frontier
V6™



FRONTIER V6 COMPACT 40

Size matters. The baby of the family, the Compact 40 is ideal for negotiating tight turns. Small but perfectly formed, the only thing that we've reduced is the size. Still the same frame, motors and seating. We just made it smaller.

"WE'VE BEEN CREATING CUSTOM MADE VEHICLES FOR 25 YEARS, SO THAT YOU CAN LIVE THE LIFE YOU CHOOSE"

Ashley Daff
Director



ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.

BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

PROUDLY AUSTRALIAN DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

Basic specifications:**

- Length: 910mm
- Width: 655mm
- Seat to floor height: 435mm
- Turning radius: 490mm
- Weight limit 182kg (400lbs), 155kg (340lbs) with lift or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Hybrid 14"x3" grey and black
- Castors: 200x50 grey and black
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Sling backrest and rehab seat pan (from 12"-24" wide)
- Crash test compliant tie down points.

Rehab positioning:

- Power seat options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and Backrests: MPS, rehab backrest, adjustable seat pans and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, swingaway, power and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any custom requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- 24V LED Lighting kit
- USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, drinks holder, oxygen and vent carriers
- Special controls: Attendant, head, chin and environmental controls all available
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:*

- Every day off-road recreation
- Urban outdoor living and working
- Roomy indoor spaces
- Tight indoor spaces.

Base options:



V6 COMPACT 40 - GREY AND BLACK WHEELS



V6 COMPACT 73 - GREY AND BLACK WHEELS



V6 HYBRID - GREY AND BLACK WHEELS

*Please remember that your wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea
**Specifications subject to change. Check online for current version
**Various model options excluded



Magic Mobility +61 3 8791 5600
sales@magicmobility.com.au

Explore with us magicmobility.com.au



Scan for agent locations