

Viewel

THE MAGIC MOBILITY JEWEL
TILT IN SPACE



Magic™
MOBILITY



BASE UNIT SPECIFICATION

- LENGTH: 1010mm
- WIDTH: 575mm
- SEAT TO FLOOR HEIGHT: 560mm (front) 460mm (rear)
- HEIGHT: 1050mm
- MAX USER WEIGHT: 100Kgs
- SEATING: MPS seat, back and headrest (460 – 610mm wide)
- TILT IN SPACE
- SWINGAWAY LEGRESTS
- TIE DOWNS

POSITIONING

- SEATING: Rehab (340 - 610mm wide), DG, custom bolsters, foam in place
- HEADRESTS: Rehab, contoured, L-type, Y-type, soft cushions
- LEGRESTS: Gas elevating
- ARMPADS: Waterfall, gutter, gel
- STRAPS AND HARNESSSES: Chest, ankle, harness, posture belts
- PADS AND SUPPORTS: knee, abduction, thoracic

OPTIONS AND ACCESSORIES

- TILT: Forwards
- FRAME: RWD, FWD, splayed castors, custom
- WHEELS: 24", 20" hub brakes, solid or pneumatic
- CASTORS: Heavy duty
- ACCESSORIES: Bag hooks, pockets, sheepskin, tables, cup holder, ventilator carrier
- POWER: Tilt, recline, legrests
- ATTENDANT: Power assist units, height adjustable handles
- INTEGRATION: aftermarket seating and power wheels
- CUSTOM ACCESSORIES: On request



AUSTRALIAN

Designed
Manufactured
Owned

AUSTRALIAN

Standards Tested

CONTACTS

Magic Mobility Pty Ltd
2/112 Browns Road
Noble Park VIC 3174 Australia

TEL +61 3 8791 5600
FAX +61 3 9701 0144
EMAIL sales@magicmobility.com.au
WEB www.magicmobility.com.au

Magic™
MOBILITY