

"FINALLY A WHEELCHAIR THAT GETS ME TO THE TOP OF MY DRIVEWAY. SOMETIMES IT'S THE LITTLE THINGS LIKE BEING ABLE TO COLLECT MY OWN MAIL THAT MAKE A DIFFERENCE"

Sue
United Kingdom



MAGIC MOBILITY FRONTIER V4 FRONT WHEEL DRIVE

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.



UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be custom designed to a fixed position. It's whatever you need.

TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.



AUSTRALIAN
MADE
-
BUILT TO LAST
-
CUSTOMISED
FOR YOUR
LIFESTYLE

MAGIC MOBILITY FRONTIER V4 FRONT WHEEL DRIVE



The Magic Mobility Front Wheel Drive* (FWD) is a great choice for those who want a more compact footprint, yet maintain excellent off-road capability. The FWD has excellent climbing ability and is highly maneuverable. Its additional benefits are wider footrest options and the ability to achieve a 90° knee angle. You can choose Off-road or Hybrid wheels to ensure you get the best wheelchair for your lifestyle.

Basic specifications:** (AT wheels and castors):

- Length: 995mm
- Width: 710mm (670mm with Hybrid wheels)
- Seat to floor height: 435mm
- Turning radius: 520mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift and tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Off-road 14" black, Hybrid 14"x3" grey and black
- Castors: Off-road and Hybrid (2.50-4)
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Sling backrest and rehab seat pan (from 12"-24" wide)
- Crash test compliant tie down points.

Rehab positioning:

- Power seat options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and Backrests: MPS, rehab backrest, adjustable seat pans and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, swingaway, power and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any custom requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (off-road wheels only)
- 24V LED Lighting kit
- 12V Accessory socket and 2 x USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, tables, drinks holder, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:

- Adventurous off-road recreation
- Country living and working
- Urban outdoor living and working
- Roomy indoor spaces.

Base options:



V4 FRONT WHEEL DRIVE - BLACK RIMS



V4 FRONT WHEEL DRIVE - HYBRID GREY WHEELS



V4 FRONT WHEEL DRIVE - HYBRID BLACK WHEELS

"WE'VE BEEN
CREATING CUSTOM
MADE VEHICLES
FOR 22 YEARS,
SO THAT YOU CAN
LIVE THE LIFE
YOU CHOOSE"

Ashley Daff
Director




ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.

BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

PROUDLY AUSTRALIAN OWNED, DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

*Please remember that your wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea
**Specifications subject to change. Check online for current version
**Various model options excluded



Visit us online to see how we change lives every day:
magicmobility.com.au

Postal Address

Magic Mobility Pty Ltd
2/112 Browns Road
Noble Park VIC 3174 Australia

Tel +61 3 8791 5600

Fax +61 3 9701 0144

Email sales@magicmobility.com.au



Join the Magic Mobility social conversations!

