

**"FINALLY A WHEELCHAIR THAT GETS ME TO THE TOP OF MY DRIVEWAY. SOMETIMES IT'S THE LITTLE THINGS LIKE BEING ABLE TO COLLECT MY OWN MAIL THAT MAKE A DIFFERENCE"**

**Sue**  
United Kingdom



# MAGIC MOBILITY FRONTIER V4 REAR WHEEL DRIVE

## FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.



## UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



## REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be custom designed to a fixed position. It's whatever you need.

## TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.

AUSTRALIAN  
MADE  
-  
BUILT TO LAST  
-  
CUSTOMISED  
FOR YOUR  
LIFESTYLE



# MAGIC MOBILITY FRONTIER V4 REAR WHEEL DRIVE



The Magic Mobility Rear Wheel Drive\* (RWD) is a small, yet high performing base. With the weight positioned over the rear wheels it delivers excellent control and tracking, especially downhill. It's a great wheelchair on the street or off the beaten path. You can choose Off-road or Hybrid wheels to ensure you get the best wheelchair for your lifestyle.

## Basic specifications:\*\* (AT wheels and castors):

- Length: 930mm
- Width: 710mm (670mm with Hybrid wheels)
- Seat to floor height: 435mm
- Turning radius: 520mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Off-road 14" black, Hybrid 14"x3" black and grey
- Castors: Off-road and Hybrid (2.50-4)
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Sling backrest and rehab seat pan (from 12"-24" wide)
- Crash test compliant tie down points.

## Rehab positioning:

- Power options: 300mm seat lift, tilt (45 degrees), power recline, lay down, power elevating legrests
- Seating and Backrests: MPS, rehab backrest, adjustable seat pans and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, swingaway, power and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any custom requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

## Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (off-road wheels only)
- 24V LED Lighting kit
- 12V Accessory socket and 2 x USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, tables, drinks holder, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

## The right wheelchair for:

- Adventurous off-road recreation
- Country living and working
- Urban outdoor living and working
- Roomy indoor spaces.

## Base options:



V4 REAR WHEEL DRIVE - BLACK RIMS



V4 REAR WHEEL DRIVE - HYBRID GREY WHEELS



V4 REAR WHEEL DRIVE - HYBRID BLACK WHEELS

"WE'VE BEEN  
CREATING CUSTOM  
MADE VEHICLES  
FOR 22 YEARS,  
SO THAT YOU CAN  
LIVE THE LIFE  
YOU CHOOSE"

**Ashley Daff**  
Director




## ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.

## BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



## ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

## PROUDLY AUSTRALIAN OWNED, DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

\*Please remember that your wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea  
\*\*Specifications subject to change. Check online for current version  
\*\*Various model options excluded



Visit us online to see how we change lives every day:  
[magicmobility.com.au](http://magicmobility.com.au)

## Postal Address

Magic Mobility Pty Ltd  
2/112 Browns Road  
Noble Park VIC 3174 Australia

**Tel** +61 3 8791 5600

**Fax** +61 3 9701 0144

**Email** [sales@magicmobility.com.au](mailto:sales@magicmobility.com.au)



Join the Magic Mobility social conversations!

