



"THIS CHAIR HAS MET AND EXCEEDED MY EXPECTATIONS. I'VE GONE MANY MILES AND HOURS ... IN THE ALMOST FOUR YEARS I'VE HAD IT. IT HAS NEVER BROKEN DOWN ON ME OR LET ME DOWN. IT'S ALWAYS GOT ME HOME, THAT'S IMPORTANT. I HAVE CONFIDENCE WHEN I LEAVE ON MY DAILY EXCURSIONS."

Gary Hays - Texas, USA

# MAGIC MOBILITY FRONTIER V6 ALL TERRAIN

AUSTRALIAN  
MADE  
-  
BUILT TO LAST  
-  
CUSTOMISED  
FOR YOUR  
LIFESTYLE



## FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.



## UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



## REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be custom designed to a fixed position. It's whatever you need.

## TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.

**Magic**  
MOBILITY™  
BECAUSE THE WORLD'S NOT FLAT

Frontier  
**V6**™

  
AUSTRALIAN  
MADE

# MAGIC MOBILITY FRONTIER V6 ALL TERRAIN



For ultimate indoor-outdoor freedom. This is it, our most popular wheelchair worldwide due to its versatility. An engineering beauty that adapts to any environment. The Frontier V6 All Terrain\* is the ideal outdoor wheelchair, yet its size makes it a capable indoor performer. Its suspension and outdoor performance lets you tackle the bumpy terrain life throws at you.

#### Basic specifications:\*\*

- Length: 1050mm
- Width: 710mm
- Seat to floor height: 435mm
- Turning radius: 580mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift and tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Off-road 14" (145/70-6) black
- Castors: 210x60mm (2.50-4) black
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Sling backrest and rehab seat pan (from 12"-24" wide)
- WC19 crash test compliant tie down points.

#### Rehab positioning:

- PPower seat options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and Backrests: MPS, rehab backrest, adjustable seat pans and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, swingaway, power and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any custom requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

#### Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (off-road wheels only)
- 24V LED Lighting kit
- USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, drinks holder, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
- Hybrid wheel kit: narrower, grey and black drive wheels and castors, when off-road capabilities are not needed. This also reduces the size of the chair further
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

#### The right wheelchair for:

- Adventurous off-road recreation
- Everyday off-road, country living and working
- Urban outdoor living and working
- Roomy indoor spaces.

#### Base options:



V6 ALL TERRAIN - BLACK RIMS



V6 HYBRID - BLACK WHEELS



V6 HYBRID - GREY WHEELS

"WE'VE BEEN  
CREATING CUSTOM  
MADE VEHICLES  
FOR 22 YEARS,  
SO THAT YOU CAN  
LIVE THE LIFE  
YOU CHOOSE"

Ashley Daff  
Director




#### ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.

#### BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



#### ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

#### PROUDLY AUSTRALIAN OWNED, DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

\*Please remember that your wheelchair has electric motors and must

NEVER be driven through water, into rivers, creeks and the sea

\*\*Specifications subject to change. Check online for current version

\*\*Various model options excluded

**Magic**  
MOBILITY  
BECAUSE THE WORLD'S NOT FLAT

Visit us online to see how  
we change lives every day:  
[magicmobility.com.au](http://magicmobility.com.au)

#### Postal Address

Magic Mobility Pty Ltd  
2/112 Browns Road  
Noble Park VIC 3174 Australia

Tel +61 3 8791 5600

Fax +61 3 9701 0144

Email [sales@magicmobility.com.au](mailto:sales@magicmobility.com.au)



Join the Magic Mobility  
social conversations!

