



**"I'M INDEPENDENT AGAIN AND ABLE TO KEEP UP WITH PEOPLE WHO AREN'T IN A WHEELCHAIR"**

**Nicole**  
Melbourne, Australia

**FULL LIGHTING PACKAGES**

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.

**UP CLOSE AND PERSONAL**

Side transfers aren't the only reason you need to flip your armrest out of the way.



**POWER OPTIONS**

Consider a seat lift to help you reach, transfer and face life head on, or power tilt for improved maneuvering over obstacles.



**TAKE CONTROL**

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.

# MAGIC MOBILITY FRONTIER V6 COMPACT C73 URBAN

AUSTRALIAN  
MADE  
-  
BUILT TO LAST  
-  
CUSTOMIZED  
FOR YOUR  
LIFESTYLE



201 Growth Parkway  
Angola, IN 46703  
United States

Tel 1-260-665-2769

Email [iim@mobility-usa.com](mailto:iim@mobility-usa.com)

Web [mobility-usa.com](http://mobility-usa.com)

V2/USA/OCT17





# MAGIC MOBILITY FRONTIER V6 URBAN C73



Live your urban life...plus the freedom to go further!  
Offering the same reliable engineering as its big brother, Frontier V6 AT, the Compact C73 Urban\* delivers best-in-class performance so you can tackle everything your day brings. Whether navigating curbs, steep hills, mulch or gravel terrains, its superior mobility and strong traction make it the perfect all-rounder, so you can get on with life and do little bit more than you dare to imagine.

#### Basic specifications:\*\*

- Length: 39.3"
- Width: 25.5"
- Seat to floor height: 17"
- Turning radius: 21"
- Weight limit: 350lbs
- Top speed: 5.5mph
- Motors: 650W Gear in Line
- Drive wheels designation: 12.5" x 2.25" or 14-3 (14" x 3")
- Front castors designation: 200 x 50
- Rear castors designation: 2.80/2.50-4 (8")
- Batteries: Group 34, 60 Ah Gel (Group 24, 75 Ah LA/AGM and Gel options)
- VR2 controller
- Center post legrest with one piece footplate
- Crash test compliant tie down points.

#### Rehab positioning:

- Power options (limited to 2 options with VR2 controllers): 12" seat lift, tilt (45 degrees), power recline, lay down, power elevating legrests, transfer tilt (8 degrees and up to 2¾" floor to seat height reduction), transfer tilt with tilt only
- Seating and backrests: MPS, rehab backrest and adjustable seat pan and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: Center post, swingaway, power elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable
- A full range of armpads, straps and supports are all available from our extensive catalogue.

#### Upgrades, options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Attendant control module
- 24V LED Lighting module kit
- 12V Accessory socket and 2 x USB charger
- Luggage rack, bags, tables, drinks holder, oxy and vent carriers
- Retractable docking pin (with R-Net upgrade only)
- Special controls: Head, chin, foot, sip and puff, Bluetooth environmental controls all available with upgrade to R-Net
- Upgradable to All-Terrain capability - see your dealer for details
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

#### The right wheelchair for:

- Urban outdoor living and working
- Roomy indoor spaces
- Tight indoor spaces.

#### Power options:



BASE MODEL



WITH POWER LEGS AND RECLINE



WITH LIFT AND TILT

"WE'VE BEEN  
CREATING CUSTOM  
MADE VEHICLES  
FOR 22 YEARS,  
SO THAT YOU CAN  
LIVE THE LIFE  
YOU CHOOSE"

Ashley Daff  
Director




#### ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.

#### BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



#### ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



#### PROUDLY AUSTRALIAN OWNED AND DESIGNED. ASSEMBLED IN THE USA

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.



Visit us online to see how we change lives every day:  
[magicmobility.com.au](http://magicmobility.com.au)

**Postal Address**  
Magic Mobility Pty Ltd  
2/112 Browns Road  
Noble Park VIC 3174 Australia

**Tel** +61 3 8791 5600  
**Fax** +61 3 9701 0144  
**Email** sales@magicmobility.com.au



Join the Magic Mobility social conversations!



\*Please remember that your wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea

\*\*Specifications subject to change. Check online for current version

\*\*Various model options excluded