

"THANK YOU FOR MAKING MY NEW CAMERA MOUNT FOR MY WHEELCHAIR. WORDS CANNOT EXPRESS THE ENORMOUS DIFFERENCE IT HAS MADE IN MY LIFE AND WORK."

Nicole  
Melbourne, VIC



# MAGIC MOBILITY FRONTIER V6 HYBRID

## FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.



## UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



## REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be custom designed to a fixed position. It's whatever you need.

## TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.



AUSTRALIAN  
MADE  
-  
BUILT TO LAST  
-  
CUSTOMISED  
FOR YOUR  
LIFESTYLE

**Magic**  
MOBILITY™  
BECAUSE THE WORLD'S NOT FLAT

Frontier  
**V6**™

  
AUSTRALIAN  
MADE

# MAGIC MOBILITY FRONTIER V6 HYBRID



Two wheelchairs in one! The Frontier V6 Hybrid is the narrower twin of the V6 All Terrain. It's an ideal indoor wheelchair with unrivaled off-road capability. For those seeking a true all terrain experience, you can also select the off-road wheel kit. Choose the best of both worlds.

## Basic specifications:

- Length: 1050mm
- Width: 695mm
- Seat to floor height: 435mm
- Turning radius: 580mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Hybrid 14" x 3" black and grey
- Castors: 210x60mm (2.50-4) grey and black
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Sling backrest and rehab seat pan (from 12-24" wide)\*
- Crash test compliant tie down points.

## Rehab positioning:

- Power options: 300mm seat lift, tilt (45 degrees), power recline, lie down, power elevating legrests, transfer tilt (8 degrees and up to 70mm floor to seat height reduction)
- Seating and Backrests: MPS, rehab backrest, adjustable seat pans and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, swingaway, power and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any custom requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

## Options and accessories:

- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (off-road wheels only)
- 24V LED Lighting kit
- USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, drinks holder, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
- Off-road wheel kit: 145/70-6 black off-road wheels, black castors, silver or black rims and jack.

## The right wheelchair for:

- Everyday off-road recreation
- Urban outdoor living and working
- Roomy indoor spaces.

## Base options:



V6 HYBRID - GREY WHEELS



V6 ALL TERRAIN - OFF-ROAD WHEELS  
BLACK RIMS



V6 ALL TERRAIN - OFF-ROAD WHEELS  
SILVER RIMS

"WE'VE BEEN  
CREATING CUSTOM  
MADE VEHICLES  
FOR 22 YEARS,  
SO THAT YOU CAN  
LIVE THE LIFE  
YOU CHOOSE"

Ashley Daff  
Director




## ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.

## BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



## ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

## PROUDLY AUSTRALIAN OWNED, DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

\*Specifications subject to change.  
Check online for current version

\*Various model options excluded

**Magic**  
MOBILITY  
BECAUSE THE WORLD'S NOT FLAT

Visit us online to see how  
we change lives every day:  
[magicmobility.com.au](http://magicmobility.com.au)

## Postal Address

Magic Mobility Pty Ltd  
2/112 Browns Road  
Noble Park VIC 3174 Australia

**Tel** +61 3 8791 5600

**Fax** +61 3 9701 0144

**Email** [sales@magicmobility.com.au](mailto:sales@magicmobility.com.au)



Join the Magic Mobility  
social conversations!

