



"WITHIN 1 HOUR OF DELIVERY I WAS IN THE BUSH AND ON THE BEACH! A DREAM COME TRUE AND PROBABLY ONE OF THE BEST FEELINGS I'VE EVER HAD, ABSOLUTE FREEDOM."

Nick
New South Wales

HIT THE ROAD

Or the sand, or the bush, or the snow.*

*Not to be driven through water

MAKING YOUR CHAIR WORK FOR YOU

Magic Mobility wheelchairs can be found on building sites, farms, stables, national parks - all allowing people to continue to work.



DO YOU KNOW YOUR LIMITS?

The Extreme X8 will help you find them, challenge yourself and embrace your freedom.



MAGIC MOBILITY EXTREME X8



**AUSTRALIAN
MADE**
—
BUILT TO LAST
—
**CUSTOMISED
FOR YOUR
LIFESTYLE**

V2/AUST/OCT17

Magic
MOBILITY™
BECAUSE THE WORLD'S NOT FLAT

Extreme
X8™



MAGIC MOBILITY EXTREME X8



The Extreme X8* is Magic Mobility's flagship powerchair. Unrivalled off-road performance and reliability make it perfect for outdoor adventurers, explorers and rural workers. The Extreme X8 is the most compact 4WD powerchair on the market.

Basic specifications:**

- Length: 1030mm
- Width: 710mm
- Seat to floor height: 470mm (490mm with swing away leg rests)
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift or tilt
- Top speed: 10kph (6.2mph)
- Motors: 700W Gear in Line
- Wheels: Off-road 14" (145/70-6) black
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Sling backrest and rehab seat pan (from 12"-24" wide).

Rehab positioning:

- Power options: 300mm seat lift, tilt (45 degrees), power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: centre post, swingaway, power and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any custom requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Puncture resistant tyre liners
- 24V LED Lighting kit
- USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, drinks holder, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
- Crash test compliant tie downs
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:

- Extreme off-road adventures
- Everyday off-road country living, working and rural productivity
- Active outer urban living.

Options:



RED - BLACK RIMS



SKY BLUE - BLACK RIMS



TANGERINE - BLACK RIMS

"WE'VE BEEN
CREATING CUSTOM
MADE VEHICLES
FOR 22 YEARS,
SO THAT YOU CAN
LIVE THE LIFE
YOU CHOOSE"

Ashley Daff
Director




ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.

BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

PROUDLY AUSTRALIAN OWNED, DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

*Please remember that your wheelchair has electric motors and must

NEVER be driven through water, into rivers, creeks and the sea

**Specifications subject to change. Check online for current version

**Various model options excluded

Magic
MOBILITY
BECAUSE THE WORLD'S NOT FLAT

Visit us online to see how
we change lives every day:
magicmobility.com.au

Postal Address

Magic Mobility Pty Ltd
2/112 Browns Road
Noble Park VIC 3174 Australia

Tel +61 3 8791 5600

Fax +61 3 9701 0144

Email sales@magicmobility.com.au



Join the Magic Mobility
social conversations!

