

"I CAN DO THINGS I NEVER THOUGHT I WOULD BE ABLE TO DO AGAIN"

Nicole
Melbourne, VIC




MAGIC MOBILITY FRONTIER V6 COMPACT 73

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.

UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.

PUT IT IN THE BOOT AND HIT THE ROAD

This folding backrest is a great feature for various backrest options. A simple quick release pin allows you to fold your wheelchair to nearly half the size.



**AUSTRALIAN
MADE**
-
BUILT TO LAST
-
**CUSTOMISED
FOR YOUR
LIFESTYLE**

Magic
MOBILITY
BECAUSE THE WORLD'S NOT FLAT

Frontier
V6

**AUSTRALIAN
MADE**

MAGIC MOBILITY FRONTIER V6 COMPACT 73



For some, the answer lies in the middle. Not the biggest, not the smallest but just right. The Frontier family has proven itself for over a decade. With the mid-sized castor arms, the Compact 73 handles curbs and takes potholes in its stride, giving you improved freedom in and around the streets.

Basic specifications:

- Length: 965mm
- Width: 655mm
- Seat to floor height: 435mm
- Turning radius: 530mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Hybrid 14x3 grey and black
- Castors: 200x50 grey and black
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Sling backrest and rehab seat pan (from 12-24" wide)*
- Crash test compliant tie down points.

Rehab positioning:

- Power options: 300mm seat lift, tilt (45 degrees), power recline, lie down, power elevating legrests, transfer tilt (8 degrees and up to 70mm floor to seat height reduction)
- Seating and Backrests: MPS, rehab backrest, adjustable seat pans and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, swingaway, power and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any custom requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

Options and accessories:

- Fold forward back
- Retractable docking pin
- 24V LED Lighting kit
- USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, drinks holder, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available.

The right wheelchair for:

- Every day off-road recreation
- Urban outdoor living and working
- Roomy indoor spaces
- Tight indoor spaces.

Base options:



V6 COMPACT 73 - GREY AND BLACK WHEELS



V6 COMPACT 40 - GREY AND BLACK WHEELS



V6 HYBRID - GREY AND BLACK WHEELS

"WE'VE BEEN CREATING CUSTOM MADE VEHICLES FOR 22 YEARS, SO THAT YOU CAN LIVE THE LIFE YOU CHOOSE"

Ashley Daff
Director




ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.

BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

PROUDLY AUSTRALIAN OWNED, DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

*Specifications subject to change. Check online for current version

*Various model options excluded

Magic
MOBILITY
BECAUSE THE WORLD'S NOT FLAT

Visit us online to see how we change lives every day:
magicmobility.com.au

Postal Address

Magic Mobility Pty Ltd
2/112 Browns Road
Noble Park VIC 3174 Australia

Tel +61 3 8791 5600

Fax +61 3 9701 0144

Email sales@magicmobility.com.au



Join the Magic Mobility social conversations!

