

FRONTIER DRIVE WHEEL OPTIONS

April 2022

FRONTIER DRIVE WHEELS



3.00-8 (14X3) SOLID OR PNEUMATIC

Polyurethane solid inserts cannot be punctured but give a harder ride characteristic

Pneumatic gives a softer ride with better shock absorption and traction

Grey or black tyres

Neverflat optional



5.30/4.50-6 CROSSOVER PNEUMATIC

Great shock absorption and traction with a less aggressive tread. Cannot be run at low pressure

Grey or black tyres

Neverflat optional



145/70-6 ALL TERRAIN PNEUMATIC

The ultimate off-road tyre, run at low pressure for great traction

Black or silver rims

Tyre liner and neverflat available

Contact Magic Mobility if you require further information

3 International Court, Scoresby Vic 3179

Tel +61 3 8791 5600

Email sales@magicmobility.com.au

<http://www.magicmobility.com.au>