

**magic**  
*Mobility*™

Magic 360

**GO YOUR  
OWN WAY**



BUILT  
TO LAST

CUSTOMISED FOR  
YOUR LIFESTYLE



AUSTRALIAN  
MADE

**360** magic™



# URBAN MANOEVRABILITY NOW GOES OFF-ROAD

## 610-660mm wide compact body

The first chair that's compact enough to move easily indoors, and equipped with everything you need for the great outdoors.

## NO JIGGLE

### Improved all-wheel suspension

Engineered by the all terrain experts, and with fully independent all-wheel suspension to take care of you and your body.



CONFIGURE YOUR CHAIR TO SUIT WHERE YOU NEED TO GO AND

WHAT YOU WANT TO DO

# 360° Freedom

The Magic 360 is the world's most versatile compact chair with off-road capacity. Tailored for the world you live in, the Magic 360 lets you choose the life you want. Configure the body style of the chair to suit your tastes. Select the wheels for the adventures you choose. Tailor the world to suit you. Your chair. Your choice.

Because the world's not flat™

FROM TIGHT INDOORS, URBAN STREETS TO OFF-ROAD



## KEEP UP THE PACE

### 5kph/3mph with power lift

There's no slowing down with 12" lift because the world doesn't stop for anyone.

## GO ANYWHERE VERSATILITY

### Choose from three interchangeable wheel options

Go where you choose, from tight indoors and offices to parks, gardens and off-road terrains.



**Urban 14" drive wheels**  
for ultimate urban manoeuvrability (610mm / 24" wide)



**Crossover 14" drive wheels**  
for added suspension and comfort (635mm / 25" wide)



**Off-road 14" drive wheels**  
for ultimate traction outdoors (660mm / 26" wide)

# 2X THE POWER TO PERFORM

## 800 watt motors

Only the most powerful motors will enable you to take on more in life with absolute confidence.

# CUSTOMISE TO YOUR LIFESTYLE

This is your life and your journey. Choose how you want it to be.

## CHOOSE YOUR > DRIVE WHEEL STYLE

Flexibility is key - and three different wheel styles open up your world like never before. Easy to change, and all are robust performers.



INTERCHANGEABLE WHEELS	URBAN 14"	CROSSOVER 14"	OFF-ROAD 14"
WIDTH OF BASE	610mm / 24"	635mm / 25"	660mm / 26"
TERRAIN	For tight indoors and life around town	For everything in between	For going further off-road

## CHOOSE YOUR > COLOUR STYLE


Take charge of your image, and choose the body colour that suits you best. From textured to high gloss. Cool elegance or vibrant shimmer to rugged tones.



## MAGIC 360 SPECIFICATIONS

USER WEIGHT LIMIT	160kg / 350lb		
MAXIMUM SPEED	10kph / 6mph		
	with 12" power lift	5kph / 3mph	
LENGTH	983mm / 38.7"		
WIDTH OF BASE	610mm / 24"	635mm / 25"	660mm / 26"
SEAT TO FLOOR HEIGHT	435mm / 17"		
GROUND CLEARANCE	90mm / 3.5"		
CASTORS (front & back)	Ø200mm / 8" solid		
TURNING RADIUS	555mm / 21"		
MOTORS	2 x 800W gear in line		
CHAIR WEIGHT (minimum)	base, batteries, rehab seating, centre post legrest, no power lift functions	140kg / 309lb	
(standard)	base, batteries, rehab seating with power lift, tilt, recline and legrest	175kg / 386lb	
DRIVE ELECTRONICS	R-net 120A		
BATTERY SIZE	73Ah gel		
SEATING (depth & width)	305-560mm / 12-22"		

## STANDARD POWER FUNCTIONS



**magic**  
Mobility™  
Because the world's not flat

Magic Mobility +61 3 8791 5600  
sales@magicmobility.com.au  
Explore with us [magicmobility.com.au](http://magicmobility.com.au)



Scan for agent locations

Please remember that your wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea. Check online for current version. The information and specifications contained in this brochure are correct and up-to-date at the time of printing and may be subject to change. ©Magic Mobility. Product of Australia.

