

TAILORED FOR YOUR LIFESTYLE

This is your life and your journey. Choose how you want it to be.



CHOOSE YOUR
> DRIVE WHEEL STYLE

Flexibility is key – three different wheel styles open up your world like never before. Easily interchangeable to suit your daily terrain.



WHEEL STYLE	URBAN 14"	CROSSOVER 14"	OFF-ROAD 14"
WIDTH OF BASE	610mm / 24"	635mm / 25"	660mm / 26"
TERRAIN	For tight indoors and life around town	For everything in between	For going further off-road

CHOOSE YOUR
> COLOUR STYLE

Your chair should reflect your style and personality. Choose from a solid metallic accent colour or a hydro-dipped pattern.



MAGIC 360
SPECIFICATIONS

USER WEIGHT LIMIT	160kg / 350lb		
MAX. SPEED	10kph / 6mph		
MAX. SPEED (POWER LIFT)	5kph / 3mph		
LENGTH OF BASE	983mm / 38.7"		
WIDTH OF BASE	610mm / 24"	635mm / 25"	660mm / 26"
SEAT TO FLOOR HEIGHT	435mm / 17"		
GROUND CLEARANCE	90mm / 3.5"		
CASTORS <i>(front & back)</i>	Ø200mm / 8" solid		
BASE TURNING RADIUS	555mm / 21"		
MOTORS	2 x 800W gear in line		
CHAIR WEIGHT <i>(minimum)</i>	base, batteries, rehab seating, centre post legrest, no power lift functions		149kg / 328lb
<i>(standard)</i>	base, batteries, rehab seating with power lift, tilt, recline and legrest		175kg / 386lb
DRIVE ELECTRONICS	R-net		
BATTERY SIZE	73Ah gel		
SEATING <i>(width)</i>	14" - 22"		
<i>(depth)</i>	12" - 24"		

STANDARD
POWER FUNCTIONS



Magic Mobility +61 3 8791 5600
enquiries@magicmobility.com.au
Explore with us magicmobility.com.au



Scan for agent locations

Please remember that your wheelchair has electric motors and must NOT be driven through water, into rivers, creeks or the sea. Check online for current version. The information and specifications contained in this brochure are correct and up-to-date at the time of printing and may be subject to change. ©Magic Mobility. Product of Australia.



M360_ENG/AUGUST26_VC10





NO JIGGLE

Improved all-wheel suspension

Engineered by the all terrain experts, and with fully independent all-wheel suspension to take care of you and your body.



KEEP UP THE PACE

5kph/3mph with up to 12" lift

Keep moving even when elevated, because the world doesn't stop for anyone.

GO ANYWHERE VERSATILITY

Choose from three interchangeable wheel options

Go where you choose, from tight indoors and offices to parks, gardens and off-road terrains.



Urban 14" drive wheels
for ultimate manoeuvrability (610mm / 24" wide)



Crossover 14" drive wheels
for added suspension and comfort (635mm / 25" wide)



Off-road 14" drive wheels
for ultimate traction outdoors (660mm / 26" wide)

CONFIGURE YOUR CHAIR TO SUIT WHERE YOU NEED TO GO AND

WHAT YOU WANT TO DO

360° Freedom

The Magic 360 is the world's most versatile compact chair with off-road capacity. Tailored for the world you live in, the Magic 360 lets you choose the life you want. Configure the body style of the chair to suit your tastes. Select the wheels for the adventures you choose. Tailor the world to suit you. Your chair. Your terrain. Your choice.

Because the world's not flat™

FROM

TIGHT INDOORS, URBAN

STREETSCAPES TO OFF-ROAD

URBAN MANOEUVRABILITY NOW GOES OFF-ROAD

610-660mm wide compact body

The first chair that's compact enough to move easily indoors, and equipped with everything you need for the great outdoors.



2X THE POWER TO PERFORM

800 watt motors

Only the most powerful motors will enable you to take on more in life with absolute confidence.