



"WITHIN ONE HOUR OF DELIVERY I WAS IN THE BUSH AND ON THE BEACH! A DREAM COME TRUE AND PROBABLY ONE OF THE BEST FEELINGS I'VE EVER HAD, ABSOLUTE FREEDOM."

Nick
Australia

magic
Mobility®

Extreme X8

**OFF-ROAD
AND OFF
THE MAP**

HIT THE ROAD

Or the sand, or the bush, or the snow.* The Extreme X8 is your partner to almost any destination.

*Not to be driven through water

MAKING YOUR CHAIR WORK FOR YOU

The Extreme X8 can be found on building sites and farms, in stables and national parks, helping people get where they need to go.



DO YOU KNOW YOUR LIMITS?

The Extreme X8 will help challenge you to embrace your freedom. It's your life, your choice and your terrain.



CE

V8/AUST/AUGUST26



Extreme
X8®



MAGIC MOBILITY EXTREME X8

The Magic Mobility Extreme X8 provides the freedom you've always wanted. The X8 is your partner to almost any destination. Sand, snow and every terrain you can throw at it. With proven off-road performance and reliability, the X8 is perfect for outdoor adventurers, explorers and rural workers. If the great outdoors is where you want to be, the Extreme X8 is waiting to take you there.

Basic specifications:*

- Length: 1030mm
- Width: 710mm
- Seat to floor height: 470mm (490mm with swingaway leg rests)
- Weight limit: 182kg (400lbs), 160kg (350lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 700W Gear in Line
- Wheels: All Terrain 14" (145/70-6) black
- Gel batteries (70Ah)
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12"-24" wide)

Rehab positioning:

- Power options: 300mm (12") seat lift, tilt (45 degrees), lift and tilt, power recline and more
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: centre post, side mount, swingaway, power, gas or manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- An extensive range of armpads, supports, restraints and pads are available

Options and accessories:

- Battery upgrade 90Ah (for up to 20% more battery range)
- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Puncture resistant tyre liners
- 24V LED lighting kit
- 12V accessory socket and 2 x USB charger
- Wheelchair jack
- Tray tables
- Luggage rack, bag hooks, bags, keep safe storage boxes
- Special controls: attendant, head, chin, foot, sip and puff, environmental controls
- Crash test compliant tie downs
- RAM® mounting systems for mobile phones, tablets, cameras and cup holders

The right wheelchair for:**

- Extreme off-road adventures
- Everyday off-road country living, working and rural productivity
- Active outer urban living

Base options:



FLAME RED - BLACK RIMS



OCEAN BLUE - SILVER RIMS



ARMY CAMO - SILVER RIMS

"WE'VE BEEN
CREATING
HIGH-QUALITY
WHEELCHAIRS
FOR 30 YEARS,
SO THAT YOU
CAN LIVE THE
LIFE YOU
CHOOSE."

Ashley Daff
R&D Director




**ATTENTION
TO DETAIL**

We know that small touches make a big difference. Our powerchair range is configurable and modifiable to suit your individual needs.

**ROBUST
QUALITY**

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



**ENGINEERED TO
THE LAST BOLT**

We listen to and value your feedback. Our team of in-house engineers work side by side with end users and therapists. We're all about innovation; if we think we can make our powerchairs better, we will.



**AUSTRALIAN
DESIGNED
AND MADE**

We are proud that our powerchair range is designed and made in our Australian manufacturing facility.

*Specifications subject to change. Check online for current version
**Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks or the sea



Magic Mobility +61 3 8791 5600
enquiries@magicmobility.com.au
Explore with us magicmobility.com.au



Scan for agent locations