

HIT THE ROAD

Or the sand, or the bush,

*Not to be driven through water

or the snow.*

WE HAVE A SIX ACRE PROPERTY AND THE EXTREME X8 ALLOWS ME TO GET TO ALL THE UNEVEN TERRAINS AND SEE PARTS OF THE FARM THAT I HAVE NOT BEEN ABLE TO ACCESS IN FOUR YEARS! THANK YOU FOR GIVING ME MY FREEDOM AND INDEPENDENCE."

Michael Trost Maryland USA

MAKING YOUR CHAIR WORK FOR YOU

Our wheelchairs are found on building sites and farms, in stables and national parks – allowing people to continue to work.

DO YOU KNOW YOUR LIMITS?

The Extreme X8 will help challenge you to embrace your freedom.



Extreme X8 OFF-ROAD AND OFF THE MAP



TAILORED FOR YOUR LIFESTYLE



.....



MAGIC MOBILITY EXTREME X8

The Magic Mobility Extreme X8 provides the freedom you've always wanted. The X8 is your partner to almost any destination. Sand, snow and every terrain you can throw at it. With unrivalled off-road performance and reliability, the X8 is perfect for outdoor adventurers, explorers and rural workers. If the great outdoors is where you want to be, the Extreme X8 is waiting to take you there.

Basic specifications:*

- Minimum base length: 40.6"
- Minimum base width: 28"
- Seat to floor height: 18.5" (19.5" with swingaway legrests)
- Weight limit: 400lbs, 340lbs with lift and/or tilt
- Top speed: 6.2mph
- Motors: 700W Gear in Line
- Wheels: Off-road 14" (145/70-6) black
- Batteries: Group 24, 75Ah LA/AGM and Gel
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12"-24" wide).

Rehab positioning:

- Power seat options: 12" lift, tilt (45 degrees), power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: centre post, side mount, swingaway, power or gas elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are available
- 24V LED Lighting kit 2 x USB charger

Options and accessories:

(CJSM 2 or JSM joystick)

Puncture resistant tyre liners

One-click activator

Fold forward back

- Jack
- Luggage rack, bag hooks, bags, gun mounts, pockets, keep safe secret boxes, tables, oxygen and vent carriers
- · Special controls: Attendant, head, chin, foot, sip and puff, environmental controls
- Crash test compliant tie downs
- RAM[®] mounting system (tablet) and swing away tablet holder, cell phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:**

- Extreme off-road adventures
- Everyday off-road country living, working and rural productivity
- Active outer urban living.

Base options:



OCEAN BLUE - BLACK AT WHEELS & RIMS



STEALTH BLACK - BLACK AT WHEELS & RIMS



FLAME RED - BLACK AT WHEELS & RIMS

"WE'VE BEEN CREATING MADE-TO-ORDER WHEELCHAIRS FOR OVER 25 YEARS, SO THAT YOU CAN LIVE THE LIFE YOU CHOOSE"

Ashley Daff R&D Director

ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team is hands on, working side by side with customers and therapists. The entire Magic team works continuously to improve the experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



BUILT TO LAST

We never compromise on

paramount.

quality. Your powerchair is your

mobility, therefore reliability is

PROUDLY **AUSTRALIAN** DESIGNED AND ASSEMBLED IN THE USA

We love that our product is made in Australia. We love meeting our customers, hearing your ever changing needs and taking it back to the drawing board.





UNRISE Sunrise Medical (US) LLC, Fresno, CA 93727 MK-129913 Rev.A©02 2021

800.333.4000 www.SunriseMedical.com



Mobility

