

“FINALLY A WHEELCHAIR THAT GETS ME TO THE TOP OF MY DRIVEWAY. SOMETIMES IT’S THE LITTLE THINGS LIKE BEING ABLE TO COLLECT MY OWN MAIL THAT MAKE A DIFFERENCE”



Sue
United Kingdom

magic
Mobility®

Frontier V4 FWD
FOR THE ROAD LESS TRAVELLED

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on.



TAKE CONTROL

An extensive range of advanced control systems and environmental controls are available. Easy to program and use.



REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily configurable to suit your individual needs.

TAILORED FOR YOU

Select the right headrest, backrest, armrests, legrests, seating and cushioning to suit your lifestyle.



V4/AUST/MARCH23

TAILORED FOR YOUR LIFESTYLE



AUSTRALIAN MADE

Frontier
V4®



FRONTIER V4 FRONT WHEEL DRIVE

A powered wheelchair for people who know what they want. The V4 FWD offers a great solution to some of life's most tricky corners. The thinking person's chair, the V4 FWD is designed for those wanting to get around their homes with ease, using desks and tables, or driving frequently. Choose All Terrain, Crossover or Hybrid wheels to ensure you get the best powerchair for your lifestyle.

"WE'VE BEEN
CREATING
HIGH-QUALITY
WHEELCHAIRS
FOR ALMOST
30 YEARS, SO
THAT YOU CAN
LIVE THE LIFE
YOU CHOOSE."

Ashley Daff
R&D Director



Basic specifications:*

- Length: 995mm
- Width: 710mm - AT wheels
700mm - Crossover wheels
640mm - Hybrid wheels
- Seat to floor height: 435mm
- Base turning radius: 520mm
- Weight limit: 182kg (400lbs),
155kg (340lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: All Terrain 14" black,
Crossover 14" black or grey,
Hybrid 14"x3" black or grey
(solid or pneumatic)
- Castors: 220mm (2.80/2.50-4),
black or grey (solid or pneumatic)
- Gel Batteries
- Penny and Giles controller with
LCD display (R-Net)
- Centre post legrest with one piece
footplate
- Rehab seat pan (from 12"-24" wide)
- Crash test compliant tie down
points

Rehab positioning:

- Power options: 300mm (12")
seat lift, tilt (45 degrees), lift and
tilt, transfer lift and tilt only, transfer
tilt with lift and tilt, power recline
and more
- Seating and backrests: MPS, rehab
backrest, adjustable seat pans and
seating
- Headrests and mounts: MPS,
rehab
- Legrests: centre post, side mount,
swingaway, power, gas or manual
elevating
- Footplates: 1 or 2 piece fixed/flip
up/angle adjustable, any specific
requirement
- An extensive range of armpads,
supports, restraints and pads are
available

The right wheelchair for:**

- Adventurous off-road recreation
- Country living and working
- Urban outdoor living and working
- Roomy indoor spaces

Options and accessories:

- One-click activator
(CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners
(All Terrain and Crossover
wheels only)
- 24V LED lighting kit
- 12V accessory socket and
2 x USB charger
- Wheelchair jack
- Tray tables
- Luggage rack, bag hooks, bags,
pockets
- Special controls: attendant,
head, chin, foot, sip and puff,
environmental controls
- RAM® mounting system (tablet
and swing away tablet holder,
mobile phone holder, swing away
camera mount, attendant cup
holder, action camera holder and
cup holder)

Base options:



SPACE BLUE BASE
BLACK AT TYRES & RIMS



LUCENT LIME BASE
BLACK CROSSOVER TYRES



BURNT ORANGE BASE
BLACK HYBRID TYRES

ATTENTION TO DETAIL

We know that small touches make a big difference. Our powerchair range is configurable and modifiable to suit your individual needs.

ROBUST QUALITY

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



ENGINEERED TO THE LAST BOLT

We listen to and value your feedback. Our team of in-house engineers work side by side with end users and therapists. We're all about innovation; if we think we can make our powerchairs better, we will.



AUSTRALIAN DESIGNED AND MADE

We are proud that our powerchair range is designed and made in our Australian manufacturing facility.



Magic Mobility +61 3 8791 5600
enquiries@magicmobility.com.au
Explore with us magicmobility.com.au



Scan for agent locations

*Specifications subject to change. Check online for current version
**Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks or the sea