

“FINALLY A WHEELCHAIR THAT GETS ME TO THE TOP OF MY DRIVEWAY. SOMETIMES IT’S THE LITTLE THINGS LIKE BEING ABLE TO COLLECT MY OWN MAIL THAT MAKE A DIFFERENCE”

Sue  
United Kingdom

### FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on.

### TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.

### REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be made-to-order.

### UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.

magic  
Mobility™

Frontier V4 FWD  
**FOR THE ROAD LESS TRAVELLED**



V6/AUST/NOV21

TAILORED FOR  
YOUR LIFESTYLE



AUSTRALIAN  
MADE

Frontier  
**V4™**



# FRONTIER V4 FRONT WHEEL DRIVE

A powered wheelchair for people who know what they want. The V4 FWD offers a great solution to some of life's most tricky corners. The thinking person's chair, the V4 FWD is designed for those wanting to get around their homes with ease, using desks and tables, or driving frequently. Choose off-road or hybrid wheels to ensure you get the best powerchair for your lifestyle.

"WE'VE BEEN  
CREATING  
MADE-TO-ORDER  
WHEELCHAIRS  
FOR OVER 25  
YEARS, SO THAT  
YOU CAN LIVE THE  
LIFE YOU  
CHOOSE"

**Ashley Daff**  
R&D Director




## ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

## QUALITY BUILT

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



## ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team is hands on, working side by side with customers and therapists. The entire Magic team works continuously to improve your entire experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



## PROUDLY AUSTRALIAN DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing your ever changing needs and taking it back to the drawing board.

### Basic specifications\* (AT wheels and castors):

- Length: 995mm
- Width: 710mm (670mm with Hybrid wheels)
- Seat to floor height: 435mm
- Turning radius: 520mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Off-road 14" black, Hybrid 14"x3" grey and black (solid or pneumatic)
- Castors: 220mm (2.80/2.50-4), black or grey (solid or pneumatic)
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12"-24" wide)
- Crash test compliant tie down points.

### Rehab positioning:

- Power options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: Center post, swingaway, power or gas and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are available.

### Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (off-road wheels only)
- 24V LED Lighting kit
- 12V Accessory socket and 2 x USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, tables, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls
- RAM® mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

### The right wheelchair for:\*\*

- Adventurous off-road recreation
- Country living and working
- Urban outdoor living and working
- Roomy indoor spaces.

### Base options:



SPACE BLUE - BLACK AT WHEELS & RIMS



BURNT ORANGE - BLACK HYBRID WHEELS (PNEUMATIC OR SOLID)



LUCENT LIME - GREY HYBRID WHEELS (PNEUMATIC OR SOLID)

\*Specifications subject to change. Check online for current version  
\*\*Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks and the sea



Because the world's not flat



FAMILY OF BRANDS

Magic Mobility +61 3 8791 5600  
sales@magicmobility.com.au

Explore with us [magicmobility.com.au](http://magicmobility.com.au)



Scan for agent locations