



"THIS CHAIR HAS MET AND EXCEEDED MY EXPECTATIONS. I'VE GONE MANY MILES AND HOURS ... IN THE TIME I'VE HAD IT. IT HAS NEVER LET ME DOWN. IT'S ALWAYS GOT ME HOME, THAT'S IMPORTANT. I HAVE CONFIDENCE WHEN I LEAVE ON MY DAILY EXCURSIONS."

Gary Hays
Texas, USA

magic
Mobility™

Frontier V4 FWD
FOR THE ROAD LESS TRAVELLED

POWER OPTIONS

Consider a seat lift to help you reach, transfer and face life head on, or power tilt for improved manoeuvring over obstacles.

TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.



FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on.

UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



V6/USA/AUGUST21

TAILORED FOR YOUR LIFESTYLE



AUSTRALIAN DESIGNED

Frontier
V4™



FRONTIER V4 FRONT WHEEL DRIVE

A powered wheelchair for people who know what they want. The V4 FWD offers a great solution to some of life's most tricky corners. The thinking person's chair, the V4 FWD is designed for people wanting to get around their homes with ease, using desks and tables, or driving frequently. Choose off-road or hybrid wheels to ensure you get the best powerchair for your lifestyle.

Basic specifications:*

(AT wheels and castors):

- Minimum base length: 39"
- Minimum base width: 28" (26" with Hybrid wheels)
- Seat to floor height: 17"
- Turning radius: 20"
- Weight limit: 400lbs, 340lbs with lift and/or tilt
- Top speed: 6.2mph
- Motors: 800W Gear in Line
- Wheels: Off-road 14" black, Hybrid 14"x3" grey and black
- Castors: 220mm (280/250-4), black or grey (solid or pneumatic)
- Batteries: Group 24, 75 Ah LA/AGM and Gel options
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12"-24" wide)
- Crash test compliant tie down points.

Rehab positioning:

- Power seat options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: Center post, swingaway, power or gas and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are available.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (off-road wheels only)
- 24V LED Lighting kit
- 2 x USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, tables, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls
- RAM® mounting system (tablet and swing away tablet holder, cell phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:**

- Adventurous off-road recreation
- Country living and working
- Urban outdoor living and working
- Roomy indoor spaces.

Base options:



OYNX BLACK - BLACK AT WHEELS & RIMS



FLAME RED - BLACK HYBRID WHEELS (PNEUMATIC OR SOLID)



ARMY GREEN - GREY HYBRID WHEELS (PNEUMATIC OR SOLID)

"WE'VE BEEN CREATING MADE-TO-ORDER WHEELCHAIRS FOR OVER 25 YEARS, SO THAT YOU CAN LIVE THE LIFE YOU CHOOSE"

Ashley Daff
R&D Director




ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

BUILT TO LAST

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team is hands on, working side by side with customers and therapists. The entire Magic team works continuously to improve your entire experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



PROUDLY AUSTRALIAN DESIGNED AND ASSEMBLED IN THE USA

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

*Specifications subject to change. Check online for current version
**Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks and the sea



Sunrise Medical (US) LLC, Fresno, CA 93727
MK-129913 Rev.A©02 2021

800.333.4000

www.SunriseMedical.com

