

**"FINALLY A WHEELCHAIR THAT GETS ME TO THE TOP OF MY DRIVEWAY. SOMETIMES IT'S THE LITTLE THINGS LIKE BEING ABLE TO COLLECT MY OWN MAIL THAT MAKE A DIFFERENCE"**

**Sue**  
United Kingdom



**magic**  
*Mobility*™

Frontier V4 RWD  
**ADVENTURE GAINED.  
NOTHING LOST**

**FULL LIGHTING PACKAGES**

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.



**UP CLOSE AND PERSONAL**

Side transfers aren't the only reason you need to flip your armrest out of the way.



**REHAB SEATING**

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be made-to-order.

**TAKE CONTROL**

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.



V5/AUST/APRIL21

**BUILT TO LAST**

**TAILORED FOR YOUR LIFESTYLE**



**AUSTRALIAN MADE**

*Frontier*  
**V4**™





# FRONTIER V4 REAR WHEEL DRIVE

A powered wheelchair for those who know what they want. The thinking person's chair, the V4 RWD is designed for those wanting to get around the streets or off-road, as well as indoors. Choose off-road or hybrid wheels to ensure you get the best powerchair for your lifestyle.

"WE'VE BEEN  
CREATING MADE-  
TO-ORDER  
VEHICLES FOR  
OVER 25 YEARS,  
SO THAT YOU CAN  
LIVE THE LIFE  
YOU CHOOSE"

Ashley Daff  
Director




## ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

## BUILT TO LAST

We never compromise on quality. Your wheelchair is your mobility, therefore reliability is paramount. We'd much rather spend time improving our products than repairing them.



## ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



## PROUDLY AUSTRALIAN DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

### Basic specifications:\*\* (AT wheels and castors):

- Length: 930mm
- Width: 710mm (670mm with Hybrid wheels)
- Seat to floor height: 435mm
- Turning radius: 520mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Off-road 14" black, Hybrid 14"x3" black and grey
- Castors: Off-road and Hybrid (2.50-4)
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12"-24" wide)
- Crash test compliant tie down points.

### Rehab positioning:

- Power options: 300mm seat lift, tilt (45 degrees), power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, swingaway, power and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

### Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (off-road wheels only)
- 24V LED Lighting kit
- 12V Accessory socket and 2 x USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, tables, drinks holder, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

### The right wheelchair for:\*

- Adventurous off-road recreation
- Country living and working
- Urban outdoor living and working
- Roomy indoor spaces.

### Base options:



FOREST GREEN - BLACK AT WHEELS & RIMS



SUNRISE YELLOW - BLACK HYBRID WHEELS (PNEUMATIC OR SOLID)



BLAZE BLUE - GREY HYBRID WHEELS (PNEUMATIC OR SOLID)

\*Please remember that your wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea  
\*\*Specifications subject to change. Check online for current version  
\*\*Various model options excluded



Magic Mobility +61 3 8791 5600  
sales@magicmobility.com.au

Explore with us [magicmobility.com.au](http://magicmobility.com.au)



Scan for agent locations