

“FINALLY A WHEELCHAIR THAT GETS ME TO THE TOP OF MY DRIVEWAY. SOMETIMES IT’S THE LITTLE THINGS LIKE BEING ABLE TO COLLECT MY OWN MAIL THAT MAKE A DIFFERENCE”

Sue
United Kingdom

FULL LIGHTING PACKAGES

Been out having too much fun and it’s getting late? Put your lights on.

TAKE CONTROL

An extensive range of advanced control systems and environmental controls are available. Easy to program and use.

REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily configurable to suit your individual needs.

TAILORED FOR YOU

Select the right headrest, backrest, armrests, legrests, seating and cushioning to suit your lifestyle.

magic
Mobility®

Frontier V4 RWD
**ADVENTURE
GAINED,
NOTHING LOST**



CE UK
CA

V4/AUST/JULY24



Frontier
V4®



FRONTIER V4 REAR WHEEL DRIVE

A powered wheelchair for people who know what they want. The thinking person's chair, the V4 RWD is designed for people wanting to get around the streets or off-road, as well as indoors. Choose Off-road, Crossover or Urban wheels to ensure you get the best powerchair for your lifestyle.

"WE'VE BEEN
CREATING
HIGH-QUALITY
WHEELCHAIRS
FOR 30 YEARS,
SO THAT YOU
CAN LIVE THE
LIFE YOU
CHOOSE."

Ashley Daff
R&D Director




ATTENTION TO DETAIL

We know that small touches make a big difference. Our powerchair range is configurable and modifiable to suit your individual needs.

ROBUST QUALITY

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



ENGINEERED TO THE LAST BOLT

We listen to and value your feedback. Our team of in-house engineers work side by side with end users and therapists. We're all about innovation; if we think we can make our powerchairs better, we will.



AUSTRALIAN DESIGNED AND MADE

We are proud that our powerchair range is designed and made in our Australian manufacturing facility.

Basic specifications:*

- Length: 930mm
- Width: 710mm - AT wheels
700mm - Crossover wheels
640mm - Urban wheels
- Seat to floor height: 435mm
- Base turning radius: 520mm
- Weight limit: 182kg (400lbs),
160kg (350lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Off-road 14" black,
Crossover 14" black or grey, Urban
14"x3" black or grey
(solid or pneumatic)
- Castors: 220mm (2.8/2.5-4),
black or grey (solid or pneumatic)
- Gel batteries
- Penny and Giles controller
with LCD display (R-Net)
- Centre post legrest with one piece
footplate
- Rehab seat pan (from 12"-24"
wide)
- Crash test compliant tie down
points

Rehab positioning:

- Power options: 300mm (12")
seat lift, tilt (45 degrees),
lift and tilt, power recline and more
- Seating and backrests: MPS, rehab
backrest, adjustable seat pans and
seating
- Headrests and mounts: MPS,
rehab
- Legrests: centre post, side mount,
swingaway, power, gas or manual
elevating
- Footplates: 1 or 2 piece fixed/flip
up/angle adjustable, any specific
requirement
- An extensive range of armpads,
supports, restraints and pads are
available

The right wheelchair for:**

- Adventurous off-road recreation
- Country living and working
- Urban outdoor living and working
- Roomy indoor spaces

Options and accessories:

- One-click activator
(CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners
(All Terrain and Crossover
wheels only)
- 24V LED lighting kit
- 12V accessory socket and
2 x USB charger
- Wheelchair jack
- Tray tables
- Luggage rack, bag hooks, bags,
pockets
- Special controls: attendant,
head, chin, foot, sip and puff,
environmental controls
- RAM® mounting system (tablet
and swing away tablet holder,
mobile phone holder, swing away
camera mount, attendant cup
holder, action camera holder and
cup holder)

Base and tyre options:



FOREST GREEN BASE
BLACK OFF-ROAD TYRES & RIMS



BLAZE BLUE BASE
BLACK CROSSOVER TYRES



SUNRISE YELLOW BASE
BLACK URBAN TYRES

*Specifications subject to change. Check online for current version
**Please remember that your wheelchair has electric motors and must
not be driven through water, into rivers, creeks or the sea



Magic Mobility +61 3 8791 5600
enquiries@magicmobility.com.au
Explore with us magicmobility.com.au



Scan for agent locations