



"THIS CHAIR HAS MET AND EXCEEDED MY EXPECTATIONS. I'VE GONE MANY MILES AND HOURS ... IN THE ALMOST FOUR YEARS I'VE HAD IT. IT HAS NEVER BROKEN DOWN ON ME OR LET ME DOWN. IT'S ALWAYS GOT ME HOME, THAT'S IMPORTANT. I HAVE CONFIDENCE WHEN I LEAVE ON MY DAILY EXCURSIONS."

Gary Hays - Texas, USA

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.



UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be made-to-order.

TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.

magic
Mobility™

Frontier V6 All Terrain
FOR ULTIMATE
OUTDOOR
ADVENTURE



V5/AUST/APRIL21

BUILT
TO LAST

TAILORED FOR
YOUR LIFESTYLE



AUSTRALIAN
MADE

Frontier
V6™



FRONTIER V6 ALL TERRAIN

Our most popular powered wheelchair worldwide, the Magic Mobility Frontier V6 All Terrain is ideal for the outdoors, yet its size makes it a capable indoor performer. If performance and versatility is what you are after, the Frontier V6 All Terrain will impress those who prefer the track to the tarmac.

Basic specifications:**

- Length: 1050mm
- Width: 710mm
- Seat to floor height: 435mm
- Turning radius: 580mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Off-road 14" (145/70-6) black
- Castors: 210x60mm (2.50-4) black
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12"-24" wide)
- Crash test compliant tie down points.

Rehab positioning:

- Power seat options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, swingaway, power and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (off-road wheels only)
- 24V LED Lighting kit
- USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, drinks holder, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
- Hybrid wheel kit: narrower, grey and black drive wheels and castors, when off-road capabilities are not needed. This also reduces the size of the chair further
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:*

- Adventurous off-road recreation
- Everyday off-road, country living and working
- Urban outdoor living and working
- Roomy indoor spaces.

Base options:



BLAZE BLUE - BLACK AT WHEELS & RIMS



ARMY GREEN - BLACK AT WHEELS & SILVER RIMS



OYNX BLACK - GREY HYBRID WHEELS (PNEUMATIC OR SOLID)

*Please remember that your wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea

**Specifications subject to change. Check online for current version

**Various model options excluded



Magic Mobility +61 3 8791 5600
sales@magicmobility.com.au

Explore with us magicmobility.com.au



Scan for agent locations

"WE'VE BEEN
CREATING MADE-
TO-ORDER
VEHICLES FOR
OVER 25 YEARS,
SO THAT YOU CAN
LIVE THE LIFE
YOU CHOOSE"

Ashley Daff
Director



ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

BUILT TO LAST

We never compromise on quality. Your wheelchair is your mobility, therefore reliability is paramount. We'd much rather spend time improving our products than repairing them.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

PROUDLY AUSTRALIAN DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.