



"THIS CHAIR HAS MET AND EXCEEDED MY EXPECTATIONS. I'VE GONE MANY MILES AND HOURS ... IN THE TIME I'VE HAD IT. IT HAS NEVER LET ME DOWN. IT'S ALWAYS GOT ME HOME, THAT'S IMPORTANT. I HAVE CONFIDENCE WHEN I LEAVE ON MY DAILY EXCURSIONS."

Gary Hays  
USA

**FULL LIGHTING PACKAGES**

Been out having too much fun and it's getting late? Put your lights on.



**TAKE CONTROL**

An extensive range of advanced control systems and environmental controls are available. Easy to program and use.



**REHAB SEATING**

Infinite options for infinite body shapes. Our rehab seating range is easily configurable to suit your individual needs.

**TAILORED FOR YOU**

Select the right headrest, backrest, armrests, legrests, seating and cushioning to suit your lifestyle.



Frontier V6 All Terrain

**FOR ULTIMATE  
OUTDOOR  
ADVENTURE**



V6/AUST/AUGUST26





# FRONTIER V6 ALL TERRAIN

Our most popular powered wheelchair worldwide, the Magic Mobility Frontier V6 All Terrain is ideal for the outdoors, yet its size makes it a capable indoor performer. If performance and versatility is what you are after, the Frontier V6 All Terrain will impress those who prefer the track to the tarmac. Choose All Terrain, Crossover or Hybrid wheels to ensure you get the best powerchair for your lifestyle.

**Basic specifications:\***

- Length: 1050mm
- Width: 710mm - AT wheels, 700mm - Crossover wheels, 695mm - Hybrid wheels
- Seat to floor height: 435mm
- Base turning radius: 580mm
- Weight limit: 182kg (400lbs), 160kg (350lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: All Terrain 14" black. Optional: Crossover 14" black or grey, Hybrid 14"x3" black or grey (solid or pneumatic)
- Castors: 220mm (2.80/2.50-4), black or grey (solid or pneumatic)
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12"-24" wide)
- Crash test compliant tie down points

**Rehab positioning:**

- Power options: 300mm (12") seat lift, tilt (45 degrees), lift and tilt, transfer lift and tilt only, transfer tilt with lift and tilt, power recline and more
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: centre post, side mount, swingaway, power, gas or manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- An extensive range of armpads, supports, restraints and pads are available

**The right wheelchair for:\*\***

- Adventurous off-road recreation
- Everyday off-road, country living and working
- Urban outdoor living and working
- Roomy indoor spaces

**Options and accessories:**

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (All Terrain and Crossover wheels only)
- 24V LED lighting kit
- 12V accessory socket and 2 x USB charger
- Wheelchair jack
- Tray tables
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes
- Special controls: attendant, head, chin, foot, sip and puff, environmental controls
- RAM® mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder)

**Alternate wheel kits:**

- Crossover wheels: shorter tread with good off road traction; the perfect in-between solution. Black or grey options
- Hybrid wheels: narrower tyres further reduce the width of the chair; when off road capability is not required. Black or grey and solid or pneumatic options

**Base options:**



BLAZE BLUE BASE  
BLACK AT TYRES & RIMS



ULTRA VIOLET BASE  
BLACK CROSSOVER TYRES



OYNX BLACK BASE  
GREY HYBRID TYRES

“WE’VE BEEN  
CREATING  
HIGH-QUALITY  
WHEELCHAIRS  
FOR 30 YEARS,  
SO THAT YOU  
CAN LIVE THE  
LIFE YOU  
CHOOSE.”

Ashley Daff  
R&D Director




**ATTENTION  
TO DETAIL**

We know that small touches make a big difference. Our powerchair range is configurable and modifiable to suit your individual needs.

**ROBUST  
QUALITY**

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



**ENGINEERED TO  
THE LAST BOLT**

We listen to and value your feedback. Our team of in-house engineers work side by side with end users and therapists. We’re all about innovation; if we think we can make our powerchairs better, we will.

**AUSTRALIAN  
DESIGNED  
AND MADE**

We are proud that our powerchair range is designed and made in our Australian manufacturing facility.

\*Specifications subject to change. Check online for current version  
\*\*Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks or the sea



Magic Mobility +61 3 8791 5600  
enquiries@magicmobility.com.au  
Explore with us [magicmobility.com.au](http://magicmobility.com.au)



Scan for agent locations