



"THIS CHAIR HAS MET AND EXCEEDED MY EXPECTATIONS. I'VE GONE MANY MILES AND HOURS ... IN THE TIME I'VE HAD IT. IT HAS NEVER LET ME DOWN. IT'S ALWAYS GOT ME HOME, THAT'S IMPORTANT. I HAVE CONFIDENCE WHEN I LEAVE ON MY DAILY EXCURSIONS."

**Gary Hays**  
Texas, USA

**magic**  
*Mobility*™

Frontier V6 All Terrain  
**FOR ULTIMATE  
OUTDOOR  
ADVENTURE**

**FULL LIGHTING  
PACKAGES**

Been out having too much fun and it's getting late? Put your lights on.

**TAKE CONTROL**

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.



**REHAB SEATING**

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be made-to-order.

**UP CLOSE  
AND PERSONAL**

Side transfers aren't the only reason you need to flip your armrest out of the way.



V6/USA/AUGUST21

TAILORED FOR  
YOUR LIFESTYLE



AUSTRALIAN  
DESIGNED

Frontier  
**V6**™



# FRONTIER V6 ALL TERRAIN

Our most popular powered wheelchair worldwide, the Magic Mobility Frontier V6 All Terrain is ideal for the outdoors, yet its size makes it a capable indoor performer. If performance and versatility is what you are after, the Frontier V6 All Terrain will impress those who prefer the track to the tarmac.

“WE’VE BEEN CREATING MADE-TO-ORDER WHEELCHAIRS FOR OVER 25 YEARS, SO THAT YOU CAN LIVE THE LIFE YOU CHOOSE”

Ashley Daff  
R&D Director




## ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

## BUILT TO LAST

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



## ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team is hands on, working side by side with customers and therapists. The entire Magic team works continuously to improve the entire experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

## PROUDLY AUSTRALIAN DESIGNED AND ASSEMBLED IN THE USA

We love that our product is made in Australia. We love meeting our customers, hearing your ever changing needs and taking it back to the drawing board.

### Basic specifications:\*

- Minimum base length: 41.5"
- Minimum base width: 28"
- Seat to floor height: 17"
- Turning radius 23"
- Weight limit: 400lbs, 340lbs with lift and/or tilt
- Top speed: 6.2mph
- Motors: 800W Gear in Line
- Wheels: Off-road 14" (145/70-6) black
- Castors: 220mm (2.80/2.50-4), black or grey (solid or pneumatic)
- Batteries: Group 24, 75 Ah LA/AGM and Gel (Group 34, 60 Ah Gel option)
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12"-24" wide)
- Crash test compliant tie down points.

### Rehab positioning:

- Power seat options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, side mount, swingaway, power or gas elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are available.

### Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (off-road wheels only)
- 24V LED Lighting kit
- 2 x USB charger
- Jack
- Luggage rack, bag hooks, bags, gun mounts, pockets, keep safe secret boxes, tables, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls
- Hybrid wheel kit: narrower, grey and black drive wheels and castors, when off-road capabilities are not needed. This further reduces the width of the chair
- RAM® mounting system (tablet and swing away tablet holder, cell phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

### The right wheelchair for:\*\*

- Adventurous off-road recreation
- Everyday off-road, country living and working
- Urban outdoor living and working
- Roomy indoor spaces.

### Base options:



SPACE BLUE - BLACK AT WHEELS & RIMS



FLAME RED - BLACK AT WHEELS & RIMS



OYNX BLACK - BLACK HYBRID WHEELS (PNEUMATIC OR SOLID)

\*Specifications subject to change. Check online for current version  
\*\*Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks and the sea



Sunrise Medical (US) LLC, Fresno, CA 93727  
MK-129913 Rev.A©02 2021

800.333.4000

www.SunriseMedical.com

