

"I CAN DO THINGS I NEVER THOUGHT I WOULD BE ABLE TO DO AGAIN"

Nicole  
Melbourne, VIC



magic  
Mobility™

### FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.

### UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



### TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.

### PUT IT IN THE BOOT AND HIT THE ROAD

The folding backrest is a great feature for various backrest options. A simple, quick release pin allows you to fold your wheelchair to nearly half the height.

Frontier V6 Compact 73  
**THE ULTIMATE URBAN COMPACT**



V6/AUST/APRIL21

BUILT TO LAST

TAILORED FOR YOUR LIFESTYLE



AUSTRALIAN MADE

Frontier  
**V6™**



# FRONTIER V6 COMPACT 73

For some, the answer lies in the middle. Not the biggest, not the smallest, but just right. The Frontier family has proven itself for over a decade. With mid-sized castor arms, the Compact 73 handles curbs and takes potholes in its stride, giving you improved freedom in and around the streets.

“WE’VE BEEN  
CREATING MADE-  
TO-ORDER  
VEHICLES FOR  
OVER 25 YEARS,  
SO THAT YOU CAN  
LIVE THE LIFE  
YOU CHOOSE”

**Ashley Daff**  
Director




## ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

## BUILT TO LAST

We never compromise on quality. Your wheelchair is your mobility, therefore reliability is paramount. We’d much rather spend time improving our products than repairing them.



## ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



## PROUDLY AUSTRALIAN DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

### Basic specifications:\*\*

- Length: 965mm
- Width: 655mm
- Seat to floor height: 435mm
- Turning radius: 530mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Hybrid 14x3 grey and black
- Castors: 200x50 grey and black
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12–24” wide)
- Crash test compliant tie down points.

### Rehab positioning:

- Power seat options: 12” lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, swingaway, power and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

### Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- 24V LED Lighting kit
- USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, drinks holder, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available.
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

### The right wheelchair for:\*

- Every day off-road recreation
- Urban outdoor living and working
- Roomy indoor spaces
- Tight indoor spaces.

### Base options:



SUNRISE YELLOW - GREY WHEELS  
(PNEUMATIC OR SOLID)



FOREST GREEN - GREY WHEELS  
(PNEUMATIC OR SOLID)



SPACE BLUE - BLACK WHEELS  
(PNEUMATIC OR SOLID)

\*Please remember that your wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea  
\*\*Specifications subject to change. Check online for current version  
\*\*Various model options excluded



Magic Mobility +61 3 8791 5600  
sales@magicmobility.com.au

Explore with us [magicmobility.com.au](http://magicmobility.com.au)



Scan for agent locations