

"I CAN DO THINGS I NEVER THOUGHT I WOULD BE ABLE TO DO AGAIN"

Nicole
Melbourne, VIC



magic
Mobility™

PUT IT IN THE BOOT AND HIT THE ROAD

The folding backrest is a great feature for various backrest options. A simple, quick release pin allows you to fold your wheelchair to nearly half the height.

UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on.

Frontier V6 Compact 73
THE ULTIMATE URBAN COMPACT



V6/AUST/NOV21

TAILORED FOR
YOUR LIFESTYLE



AUSTRALIAN
MADE

Frontier
V6™



FRONTIER V6 COMPACT 73

For some, the answer lies in the middle. Not the biggest, not the smallest, but just right. The Frontier family has proven itself for over a decade. With mid-sized castor arms, the Compact 73 handles kerbs and takes potholes in its stride, giving you improved freedom in and around the streets.

“WE’VE BEEN
CREATING
MADE-TO-ORDER
WHEELCHAIRS
FOR OVER 25
YEARS, SO THAT
YOU CAN LIVE THE
LIFE YOU
CHOOSE”

Ashley Daff
R&D Director




ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

QUALITY BUILT

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team is hands on, working side by side with customers and therapists. The entire Magic team works continuously to improve your entire experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



PROUDLY AUSTRALIAN DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing your ever changing needs and taking it back to the drawing board.

Basic specifications:*

- Length: 965mm
- Width: 655mm
- Seat to floor height: 435mm
- Turning radius: 530mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Hybrid 14x3" grey and black (solid or pneumatic)
- Castors: 200x50mm grey or black (solid or pneumatic)
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12–24" wide)
- Crash test compliant tie down points.

Rehab positioning:

- Power options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: Center post, swingaway, power or gas and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are available.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- 24V LED Lighting kit
- 12V Accessory socket and 2 x USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls
- RAM® mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:**

- Every day off-road recreation
- Urban outdoor living and working
- Roomy indoor spaces
- Tight indoor spaces.

Base options:



SUNRISE YELLOW - GREY WHEELS
(PNEUMATIC OR SOLID)



FOREST GREEN - GREY WHEELS
(PNEUMATIC OR SOLID)



SPACE BLUE - BLACK WHEELS
(PNEUMATIC OR SOLID)

*Specifications subject to change. Check online for current version
**Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks and the sea



Because the world's not flat



FAMILY OF BRANDS

Magic Mobility +61 3 8791 5600
sales@magicmobility.com.au

Explore with us magicmobility.com.au



Scan for agent locations