

**"I'M  
 INDEPENDENT  
 AGAIN AND ABLE  
 TO KEEP UP WITH  
 PEOPLE WHO  
 AREN'T IN A  
 WHEELCHAIR"**

**Nicole**  
 Melbourne, Australia



**FULL LIGHTING  
 PACKAGES**

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.

**UP CLOSE  
 AND PERSONAL**

Side transfers aren't the only reason you need to flip your armrest out of the way.



**POWER OPTIONS**

Consider a seat lift to help you reach, transfer and face life head on, or power tilt for improved maneuvering over obstacles.

**TAKE CONTROL**

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.



Frontier V6 Compact C73 Urban

**THE ULTIMATE  
 URBAN  
 COMPACT**



**Exclusive USA distributor since 1999**

201 Growth Parkway  
 Angola, IN 46703  
 United States

**Tel** 1-260-665-2769

**Email** [iim@mobility-usa.com](mailto:iim@mobility-usa.com)

**Web** [mobility-usa.com](http://mobility-usa.com)



2016/USA/OCT20

**BUILT  
 TO LAST**

**CUSTOMISED FOR  
 YOUR LIFESTYLE**



**AUSTRALIAN  
 MADE**



**Innovation  
 In Motion**





# FRONTIER V6 COMPACT 73 URBAN

For some, the answer lies in the middle. Not the biggest, not the smallest but just right. The Frontier family has proven itself for over a decade. With the mid-sized castor arms, the Compact 73 Urban handles curbs and takes potholes in its stride, giving you improved freedom in and around the streets.

### Basic specifications:\*\*

- Length: 39.3"
- Width: 25.5"
- Seat to floor height: 17"
- Turning radius: 21"
- Weight limit: 350lbs
- Top speed: 5.5mph
- Motors: 650W Gear in Line
- Drive wheels designation: 12.5" x 2.25" or 14-3 (14" x 3")
- Front castors designation: 200 x 50
- Rear castors designation: 2.80/2.50-4 (8")
- Batteries: Group 34, 60 Ah Gel (Group 24, 75 Ah LA/AGM and Gel options)
- VR2 controller
- Center post legrest with one piece footplate
- Crash test compliant tie down points.

### Rehab positioning:

- Power options (limited to 2 options with VR2 controllers): 12" seat lift, tilt (45 degrees), power recline, lay down, power elevating legrests, transfer tilt (8 degrees and up to 2¾" floor to seat height reduction), transfer tilt with tilt only
- Seating and backrests: MPS, rehab backrest and adjustable seat pan and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: Center post, swingaway, power elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable
- A full range of armpads, straps and supports are all available from our extensive catalogue.

### Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Attendant control module
- 24V LED Lighting module kit
- 12V Accessory socket and 2 x USB charger
- Luggage rack, bags, tables, drinks holder, oxy and vent carriers
- Retractable docking pin (with R-Net upgrade only)
- Special controls: Head, chin, foot, sip and puff, Bluetooth environmental controls all available with upgrade to R-Net
- Upgradable to All-Terrain capability - see your dealer for details
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

### The right wheelchair for:\*

- Urban outdoor living and working
- Roomy indoor spaces
- Tight indoor spaces.

### Base options:



BASE MODEL



WITH POWER LEGS AND RECLINE



WITH LIFT AND TILT

"WE'VE BEEN  
CREATING CUSTOM  
MADE VEHICLES  
FOR 25 YEARS, SO  
THAT YOU CAN LIVE  
THE LIFE YOU  
CHOOSE"

Ashley Daff  
Director




### ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.



### BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



### ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



### PROUDLY AUSTRALIAN DESIGNED AND ASSEMBLED IN THE USA

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.



Magic Mobility +61 3 8791 5600  
sales@magicmobility.com.au

Explore with us [magicmobility.com.au](http://magicmobility.com.au)



Scan for agent locations