

"I CAN DO THINGS I NEVER THOUGHT I WOULD BE ABLE TO DO AGAIN."

Nicole
Australia



magic
Mobility®

Frontier V6 Compact

THE ULTIMATE
URBAN
COMPACT

PACK UP AND HIT THE ROAD

The Fold Forward Backrest is a great feature when travelling with your powerchair. A simple quick release mechanism allows you to fold your chair to nearly half the height.

TAILORED FOR YOU

Select the right headrest, backrest, armrests, legrests, seating and cushioning to suit your lifestyle.



TAKE CONTROL

An extensive range of advanced control systems and environmental controls are available. Easy to program and use.

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on.



CE UK
CA

V6/AUST/AUG24



Frontier
V6®



FRONTIER V6 COMPACT

The narrow footprint of the Frontier V6 Compact provides users with the ultimate in manoeuvrability, particularly in tight indoor spaces. The Frontier family has proven itself for two decades. The Frontier family of powerchairs handles kerbs and takes potholes in its stride, giving you improved freedom in and around the streets.

“WE’VE BEEN CREATING HIGH-QUALITY WHEELCHAIRS FOR 30 YEARS, SO THAT YOU CAN LIVE THE LIFE YOU CHOOSE.”

Ashley Daff
R&D Director




ATTENTION TO DETAIL

We know that small touches make a big difference. Our powerchair range is configurable and modifiable to suit your individual needs.

ROBUST QUALITY

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



ENGINEERED TO THE LAST BOLT

We listen to and value your feedback. Our team of in-house engineers work side by side with end users and therapists. We’re all about innovation; if we think we can make our powerchairs better, we will.



AUSTRALIAN DESIGNED AND MADE

We are proud that our powerchair range is designed and made in our Australian manufacturing facility.

Basic specifications:*

- Length: 955mm
- Width: 655mm
- Seat to floor height: 435mm
- Base turning radius: 530mm
- Weight limit: 182kg (400lbs), 160kg (350lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Urban 14x3" black or grey (solid or pneumatic)
- Castors: 200x50mm black or grey (solid or pneumatic)
- Gel batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12-24" wide)
- Crash test compliant tie down points

Rehab positioning:

- Power options: 300mm (12") seat lift, tilt (45 degrees), lift and tilt, transfer lift and tilt only, transfer tilt with lift and tilt, power recline and more
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: center post, swingaway, power or gas and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are available

The right wheelchair for:**

- Every day off-road recreation
- Urban outdoor living and working
- Roomy indoor spaces
- Tight indoor spaces

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- 24V LED lighting kit
- 12V accessory socket and 2 x USB charger
- Wheelchair jack
- Tray tables
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes
- Special controls: attendant, head, chin, foot, sip and puff, environmental controls
- RAM® mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder)

Base options:



SUNRISE YELLOW BASE
GREY URBAN TYRES



FOREST GREEN BASE
GREY URBAN TYRES



SPACE BLUE BASE
BLACK URBAN TYRES

*Specifications subject to change. Check online for current version
**Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks or the sea



Magic Mobility +61 3 8791 5600
enquiries@magicmobility.com.au
Explore with us magicmobility.com.au



Scan for agent locations