



"THIS CHAIR HAS MET AND EXCEEDED MY EXPECTATIONS. I'VE GONE MANY MILES AND HOURS ... IN THE ALMOST FOUR YEARS I'VE HAD IT. IT HAS NEVER BROKEN DOWN ON ME OR LET ME DOWN. IT'S ALWAYS GOT ME HOME, THAT'S IMPORTANT. I HAVE CONFIDENCE WHEN I LEAVE ON MY DAILY EXCURSIONS."

Gary Hays - Texas, USA

magic
Mobility™

Frontier V6 FWD
FOR THE ROAD LESS TRAVELLED

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.

UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



POWER OPTIONS

Consider a seat lift to help you reach, transfer and face life head on, or power tilt for improved maneuvering over obstacles.

TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.



 Innovation
In Motion

Exclusive USA distributor since 1999

201 Growth Parkway
Angola, IN 46703
United States

Tel 1-260-665-2769

Email iim@mobility-usa.com

Web mobility-usa.com

V3b/USA/DEC19

BUILT TO LAST

CUSTOMISED FOR YOUR LIFESTYLE



AUSTRALIAN MADE

 Innovation
In Motion

Frontier
V6™



FRONTIER V6 FRONT WHEEL DRIVE

A powered wheelchair for those who know what they want. The V6 FWD offers a great solution to some of life's most tricky corners. The thinking person's chair, the V6 FWD is designed for those who are asking themselves the question about getting around their homes with ease, using desks and tables or driving frequently. Choose off-road or hybrid wheels to ensure you get the best powerchair for your lifestyle.

Basic specifications:** (AT wheels and castors):

- Length: 39"
- Width: 28" (26" with Hybrid wheels)
- Seat to floor height: 17"
- Turning radius: 20"
- Weight limit: 400lbs, 340lbs with lift and tilt
- Top speed: 6.2mph
- Motors: 800W Gear in Line
- Wheels: Off-road 14" black, Hybrid 14"x3" grey and black
- Castors: Off-road and Hybrid 2.50-4 (8")
- Batteries: Group 24, 75 Ah LA/AGM and Gel options
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Sling backrest and rehab seat pan (from 12"-24" wide)
- Crash test compliant tie down points.

Rehab positioning:

- Power seat options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and Backrests: MPS, rehab backrest, adjustable seat pans and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: Center post, swingaway, power and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any custom requirement. Articulating footplate for center legrest (in conjunction with transfer tilt)
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (off-road wheels only)
- 24V LED Lighting kit
- 12V Accessory socket and 2 x USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, tables, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:*

- Adventurous off-road recreation
- Country living and working
- Urban outdoor living and working
- Roomy indoor spaces.

Base options:



V6 FRONT WHEEL DRIVE - HYBRID GREY WHEELS



V6 FRONT WHEEL DRIVE - HYBRID BLACK WHEELS

*Please remember that your wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea

**Specifications subject to change. Check online for current version

**Various model options excluded



Magic Mobility +61 3 8791 5600
sales@magicmobility.com.au

Explore with us magicmobility.com.au



Scan for agent locations

"WE'VE BEEN
CREATING CUSTOM
MADE VEHICLES
FOR 25 YEARS, SO
THAT YOU CAN LIVE
THE LIFE YOU
CHOOSE"

Ashley Daff
Director




ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.

BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



PROUDLY AUSTRALIAN DESIGNED AND ASSEMBLED IN THE USA

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.