

"THANK YOU FOR MAKING MY NEW CAMERA MOUNT FOR MY WHEELCHAIR. WORDS CANNOT EXPRESS THE ENORMOUS DIFFERENCE IT HAS MADE IN MY LIFE AND WORK."



Nicole
Melbourne, VIC

magic
Mobility™

Frontier V6 Hybrid

**FOR THE
BEST OF
BOTH WORLDS**

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.



UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be made-to-order.

TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.



V6/AUST/APRIL21

BUILT
TO LAST

TAILORED FOR
YOUR LIFESTYLE



AUSTRALIAN
MADE

Frontier
V6™



FRONTIER V6 HYBRID

Two wheelchairs in one! The Frontier V6 Hybrid is the narrower twin of the V6 All Terrain. It's an ideal indoor wheelchair with unrivalled off-road capability. For those seeking a true all terrain experience, you can also select the off-road wheel kit. Choose the best of both worlds.

“WE’VE BEEN CREATING MADE-TO-ORDER VEHICLES FOR OVER 25 YEARS, SO THAT YOU CAN LIVE THE LIFE YOU CHOOSE”

Ashley Daff
Director




ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

BUILT TO LAST

We never compromise on quality. Your wheelchair is your mobility, therefore reliability is paramount. We'd much rather spend time improving our products than repairing them.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



PROUDLY AUSTRALIAN DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

Basic specifications:**

- Length: 1050mm
- Width: 695mm
- Seat to floor height: 435mm
- Turning radius: 580mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Hybrid 14" x 3" black and grey
- Castors: 210x60mm (2.50-4) grey and black
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12-24" wide)
- Crash test compliant tie down points.

Rehab positioning:

- Power seat options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, swingaway, power and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (off-road wheels only)
- 24V LED Lighting kit
- USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, drinks holder, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
- Off-road wheel kit: 145/70-6 black off-road wheels, black castors, silver or black rims and jack
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:*

- Everyday off-road recreation
- Urban outdoor living and working
- Roomy indoor spaces.

Base options:



OYNX BLACK - BLACK HYBRID WHEELS (PNEUMATIC OR SOLID)



BURNT ORANGE - GREY HYBRID WHEELS (PNEUMATIC OR SOLID)



ARMY GREEN - BLACK AT WHEELS & RIMS

*Please remember that your wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea
**Specifications subject to change. Check online for current version
**Various model options excluded



Magic Mobility +61 3 8791 5600
sales@magicmobility.com.au

Explore with us magicmobility.com.au



Scan for agent locations