"THANK YOU FOR MAKING MY NEW CAMERA MOUNT FOR MY WHEELCHAIR. WORDS CANNOT EXPRESS THE ENORMOUS DIFFERENCE IT HAS MADE IN MY LIFE AND WORK."



Nicole Melbourne, <u>VIC</u>

FULL LIGHTING PACKAGES

A full range of

Been out having too much fun and it's getting late? Put your lights on. A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.

TAKE CONTROL



REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be made-to-order.

UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



Frontier V6 Hybrid FOR THE BEST OF BOTH WORLDS







AUSTRALIAN MADE



Basic specifications:*

• Seat to floor height: 435mm

Weight limit: 182kg (400lbs),

Top speed: 10kph (6.2mph)

• Motors: 800W Gear in Line

• Wheels: Hybrid 14" x 3" black

and grey (solid or pneumatic)

• Castors: 220mm (2.80/2.50-4),

• Penny and Giles controller with

• Centre post legrest with one

Crash test compliant tie down

*Specifications subject to change. Check online for current version

not be driven through water, into rivers, creeks and the sea

**Please remember that your wheelchair has electric motors and must

LCD display (R-Net)

piece footplate

• Rehab seat pan

points

(from 12-24" wide)

black or grey (solid or pneumatic)

155kg (340lbs) with lift and/or tilt

• Turning radius: 580mm

Length: 1050mm

• Width: 695mm

Gel Batteries

FRONTIER V6 HYBRID

Two wheelchairs in one! The Frontier V6 Hybrid is the narrower twin of the V6 All Terrain. It's an ideal indoor wheelchair with unrivalled off-road capability. For those seeking a true all terrain experience, you can also select the off-road wheel kit. Choose the best of both worlds.

"WE'VE BEEN CREATING MADE-TO-ORDER WHEELCHAIRS FOR OVER 25 YEARS, SO THAT YOU CAN LIVE THE LIFE YOU CHOOSE"

Ashley Daff **R&D** Director



ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

QUALITY BUILT

We never compromise on



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team is hands on, working side by side with customers and therapists. The entire Magic team works continuously to improve your entire experience, from choosing a Magic Mobility wheelchair to its performance in the places

quality. Your powerchair is your mobility, therefore reliability is paramount.



Base options:



you want to go.

PROUDLY **AUSTRALIAN** DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing your ever changing needs and taking it back to the drawing board.

Rehab positioning:

• Power options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests

- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: Center post, swingaway, power or gas and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are available.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- 24V LED Lighting kit
- 12V Accessory socket and 2 x USB charger
- Jack
- · Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls
- Off-road wheel kit: 145/70-6 black off-road wheels, black castors,
- Puncture resistant tyre liners (off-road wheels only)
- RAM[®] mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:**

- Everyday off-road recreation
- Roomy indoor spaces.



Urban outdoor living and working



Because the world's not flat



Magic Mobility +61 3 8791 5600 sales@magicmobility.com.au





ARMY GREEN - BLACK AT WHEELS & RIMS

BURNT ORANGE - GREY HYBRID WHEELS

(PNEUMATIC OR SOLID)



