

“THANK YOU FOR MAKING MY NEW CAMERA MOUNT FOR MY WHEELCHAIR. WORDS CANNOT EXPRESS THE ENORMOUS DIFFERENCE IT HAS MADE IN MY LIFE AND WORK.”

Nicole
Melbourne, VIC



magic
Mobility™

Frontier V6 Hybrid

FOR THE
BEST OF
BOTH WORLDS

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.



UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be custom designed to a fixed position. It's whatever you need.

TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.



 Innovation
In Motion

Exclusive USA distributor since 1999

201 Growth Parkway
Angola, IN 46703
United States

Tel 1-260-665-2769

Email iim@mobility-usa.com

Web mobility-usa.com

V6/USA/DEC19

BUILT TO LAST

CUSTOMISED FOR YOUR LIFESTYLE



AUSTRALIAN MADE



Innovation In Motion

Frontier
V6™



FRONTIER V6 HYBRID

Two wheelchairs in one! The Frontier V6 Hybrid is the narrower twin of the V6 All Terrain. It's an ideal indoor wheelchair with unrivaled off-road capability. For those seeking a true all terrain experience, you can also select the off-road wheel kit. Choose the best of both worlds.

Basic specifications:**

- Length: 41½"
- Width: 26½"
- Seat to floor height: 17"
- Turning radius: 23"
- Weight limit: 400lbs, 340lbs with lift or tilt
- Top speed: 6.2mph
- Motors: 800W Gear in Line
- Wheels: Hybrid 14"x3" grey and black
- Castors: 2.50-4 grey and black
- Batteries: Group 24, 75 Ah LA/AGM and Gel (Group 34, 60 Ah option)
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Sling backrest and rehab seat pan (from 12"-23" wide)
- Crash test compliant tie down points.

Rehab positioning:

- Power options: 12" seat lift, tilt (45 degrees), power recline, lay down, power elevating legrests, transfer tilt (8 degrees and up to 2¾" floor to seat height reduction), transfer tiltwith tilt only
- Seating and Backrests: MPS, rehab backrest and adjustable seat pans and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, side mount, swingaway, power elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any custom requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- 24V LED Lighting kit
- USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, drinks holder, oxy and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
- Off-road wheel kit: 145/70-6 black off-road wheels, black castors and jack
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

Base options:



V6 ALL TERRAIN - OFF-ROAD WHEELS

The right wheelchair for:*

- Everyday off-road recreation
- Urban outdoor living and working
- Roomy indoor spaces.

*Please remember that your wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea
 **Specifications subject to change. Check online for current version
 **Various model options excluded

"WE'VE BEEN CREATING CUSTOM MADE VEHICLES FOR 25 YEARS, SO THAT YOU CAN LIVE THE LIFE YOU CHOOSE"

Ashley Daff
Director




ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.

BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



PROUDLY AUSTRALIAN DESIGNED AND ASSEMBLED IN THE USA

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.



Magic Mobility +61 3 8791 5600
sales@magicmobility.com.au

Explore with us magicmobility.com.au



Scan for agent locations