

"I CAN DO THINGS I NEVER THOUGHT I WOULD BE ABLE TO DO AGAIN."

Nicole
Australia



magic
Mobility®

Frontier V6 Urban

VERSATILITY FOR YOUR TERRAIN

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on.



TAKE CONTROL

An extensive range of advanced control systems and environmental controls are available. Easy to program and use.



REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily configurable to suit your individual needs.

TAILORED FOR YOU

Select the right headrest, backrest, armrests, legrests, seating and cushioning to suit your lifestyle.



CE UK
CA

V6/AUST/AUGUST26

Frontier
V6®



FRONTIER V6 URBAN

The Frontier V6 Urban is the narrower version of the V6 All Terrain. An ideal indoor powerchair with unrivalled off-road capabilities, depending on your choice of wheels. For those seeking a true all terrain experience, select the All Terrain wheel kit, or choose the Crossover option for the perfect in-between solution.

Basic specifications:*

- Length: 1050mm
- Width: 695mm - Urban wheels
700mm - Crossover wheels
710mm - AT wheels
- Seat to floor height: 435mm
- Base turning radius: 580mm
- Weight limit: 182kg (400lbs),
160kg (350lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Urban 14" x 3" black or grey (solid or pneumatic)
Optional: Crossover 14" black or grey, All Terrain 14" black
- Castors: 220mm (2.80/2.50-4), black or grey (solid or pneumatic)
- Gel batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12-24" wide)
- Crash test compliant tie down points

Rehab positioning:

- Power options: 300mm (12") seat lift, tilt (45 degrees), lift and tilt, transfer lift and tilt only, transfer tilt with lift and tilt, power recline and more
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: centre post, side mount, swingaway, power, gas or manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- An extensive range of armpads, supports, restraints and pads are available

The right wheelchair for:**

- Everyday off-road recreation
- Urban outdoor living and working
- Roomy indoor spaces

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- 24V LED lighting kit
- 12V accessory socket and 2 x USB charger
- Wheelchair jack
- Tray tables
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes
- Special controls: attendant, head, chin, foot, sip and puff, environmental controls
- Puncture resistant tyre liners (All Terrain or Crossover wheels only)
- RAM® mounting system (tablet and swingaway tablet holder, mobile phone holder, swingaway camera mount, attendant cup holder, action camera holder and cup holder)

Alternate wheel kits:

- Crossover wheels: shorter tread with good off-road traction; the perfect in-between solution. Black or grey options
- All Terrain wheel kit: 145/70-6 black off-road wheels, black castors, silver or black rims and jack

Base options:



OYNX BLACK BASE
BLACK URBAN TYRES



ULTRA VIOLET BASE
BLACK CROSSOVER TYRES



ARMY GREEN BASE
BLACK AT TYRES & RIMS

"WE'VE BEEN
CREATING
HIGH-QUALITY
WHEELCHAIRS
FOR 30 YEARS,
SO THAT YOU
CAN LIVE THE
LIFE YOU
CHOOSE."

Ashley Daff
R&D Director




**ATTENTION
TO DETAIL**

We know that small touches make a big difference. Our powerchair range is configurable and modifiable to suit your individual needs.

**ROBUST
QUALITY**

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



**ENGINEERED TO
THE LAST BOLT**

We listen to and value your feedback. Our team of in-house engineers work side by side with end users and therapists. We're all about innovation; if we think we can make our powerchairs better, we will.

**AUSTRALIAN
DESIGNED
AND MADE**

We are proud that our powerchair range is designed and made in our Australian manufacturing facility.

*Specifications subject to change. Check online for current version
**Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks or the sea



Magic Mobility +61 3 8791 5600
enquiries@magicmobility.com.au
Explore with us magicmobility.com.au



Scan for agent locations