

**magic**  
*Mobility*<sup>®</sup>

Frontier V6 Urban

# VERSATILITY FOR YOUR TERRAIN



*Frontier*  
**V6**<sup>®</sup>

**"I CAN DO  
THINGS I  
NEVER  
THOUGHT  
I WOULD BE  
ABLE TO  
DO AGAIN."**

Nicole  
Australia



## FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on.



## TAKE CONTROL

An extensive range of advanced control systems and environmental controls are available. Easy to program and use.



## REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily configurable to suit your individual needs.

## TAILORED FOR YOU

Select the right headrest, backrest, armrests, legrests, seating and cushioning to suit your lifestyle.

CE UK  
CA

V6/AUST/JULY24





# FRONTIER V6 URBAN

The Frontier V6 Urban is the narrower version of the V6 All Terrain. An ideal indoor powerchair with unrivalled off-road capabilities, depending on your choice of wheels. For those seeking a true all terrain experience, select the All Terrain wheel kit, or choose the Crossover option for the perfect in-between solution.

“WE’VE BEEN CREATING HIGH-QUALITY WHEELCHAIRS FOR 30 YEARS, SO THAT YOU CAN LIVE THE LIFE YOU CHOOSE.”

Ashley Daff  
R&D Director




## ATTENTION TO DETAIL

We know that small touches make a big difference. Our powerchair range is configurable and modifiable to suit your individual needs.

## ROBUST QUALITY

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



## ENGINEERED TO THE LAST BOLT

We listen to and value your feedback. Our team of in-house engineers work side by side with end users and therapists. We’re all about innovation; if we think we can make our powerchairs better, we will.



## AUSTRALIAN DESIGNED AND MADE

We are proud that our powerchair range is designed and made in our Australian manufacturing facility.

### Basic specifications:\*

- Length: 1050mm
- Width: 695mm - Urban wheels  
700mm - Crossover wheels  
710mm - AT wheels
- Seat to floor height: 435mm
- Base turning radius: 580mm
- Weight limit: 182kg (400lbs),  
160kg (350lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Urban 14" x 3" black or grey (solid or pneumatic)  
Optional: Crossover 14" black or grey, All Terrain 14" black
- Castors: 220mm (2.80/2.50-4), black or grey (solid or pneumatic)
- Gel batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12-24" wide)
- Crash test compliant tie down points

### Rehab positioning:

- Power options: 300mm (12") seat lift, tilt (45 degrees), lift and tilt, transfer lift and tilt only, transfer tilt with lift and tilt, power recline and more
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: centre post, side mount, swingaway, power, gas or manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- An extensive range of armpads, supports, restraints and pads are available

### The right wheelchair for:\*\*

- Everyday off-road recreation
- Urban outdoor living and working
- Roomy indoor spaces

### Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- 24V LED lighting kit
- 12V accessory socket and 2 x USB charger
- Wheelchair jack
- Tray tables
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes
- Special controls: attendant, head, chin, foot, sip and puff, environmental controls
- Puncture resistant tyre liners (All Terrain or Crossover wheels only)
- RAM® mounting system (tablet and swingaway tablet holder, mobile phone holder, swingaway camera mount, attendant cup holder, action camera holder and cup holder)

### Alternate wheel kits:

- Crossover wheels: shorter tread with good off-road traction; the perfect in-between solution. Black or grey options
- All Terrain wheel kit: 145/70-6 black off-road wheels, black castors, silver or black rims and jack

### Base options:



OYNX BLACK BASE  
BLACK URBAN TYRES



ULTRA VIOLET BASE  
BLACK CROSSOVER TYRES



ARMY GREEN BASE  
BLACK AT TYRES & RIMS

\*Specifications subject to change. Check online for current version  
\*\*Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks or the sea



Magic Mobility +61 3 8791 5600  
enquiries@magicmobility.com.au  
Explore with us [magicmobility.com.au](http://magicmobility.com.au)



Scan for agent locations